



CAT. NO. SS	CAT. NO. TT	LCP PROXIMAL TIBIA PLATE (RAFT PLATE) 3.5 mm
346.104L	646.104L	LCP PROXIMAL TIBIA PLATE LEFT 4HOLES
346.104R	646.104R	LCP PROXIMAL TIBIA PLATE RIGHT 4HOLES
346.106L	646.106L	LCP PROXIMAL TIBIA PLATE LEFT 6HOLES
346.106R	646.106R	LCP PROXIMAL TIBIA PLATE RIGHT 6HOLES
346.108L	646.108L	LCP PROXIMAL TIBIA PLATE LEFT 8HOLES
346.108R	646.108R	LCP PROXIMAL TIBIA PLATE RIGHT 8HOLES
346.110L	646.110L	LCP PROXIMAL TIBIA PLATE LEFT 10HOLES
346.110R	646.110R	LCP PROXIMAL TIBIA PLATE RIGHT 10HOLES
346.112L	646.112L	LCP PROXIMAL TIBIA PLATE LEFT 12HOLES
346.112R	646.112R	LCP PROXIMAL TIBIA PLATE RIGHT 12HOLES
346.114L	646.114L	LCP PROXIMAL TIBIA PLATE LEFT 14HOLES
346.114R	646.114R	LCP PROXIMAL TIBIA PLATE RIGHT 14HOLES

Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plate available with 4,6,8,10,12 & 14 holes.

Indication for use

- Wise Lock 3.5mm Proximal Tibia Plate is indicated for fractures of the proximal tibia, including: simple, comminuted, lateral wedge, depression, medial wedge, bicondylar combination of lateral wedge and depression, periprosthetic, and fractures with associated shaft fractures.

The plates can also be used for the treatments of nonunions, malunions, and osteopenic bone.3.5mm Wise Lock Medial Proximal Tibial Plate