



CAT. NO.	METAL	LCP Screw 2.0 mm Star Head
626.950	TT	LC - LCP Screw 2.0 mm Star Head,L- 10 mm
626.952	TT	LC - LCP Screw 2.0 mm Star Head,L- 12 mm
626.954	TT	LC - LCP Screw 2.0 mm Star Head,L- 14 mm
626.956	TT	LC - LCP Screw 2.0 mm Star Head,L- 16 mm
626.958	TT	LC - LCP Screw 2.0 mm Star Head,L- 18 mm
626.960	TT	LC - LCP Screw 2.0 mm Star Head,L- 20 mm
626.962	TT	LC - LCP Screw 2.0 mm Star Head,L- 22 mm
626.964	TT	LC - LCP Screw 2.0 mm Star Head,L- 24 mm
626.966	TT	LC - LCP Screw 2.0 mm Star Head,L- 26 mm
626.968	TT	LC - LCP Screw 2.0 mm Star Head,L- 28 mm
626.970	TT	LC - LCP Screw 2.0 mm Star Head,L- 30 mm

Specification

- Available in Titanium.
- LCP Screw 2.0 mm Star Head Length ranges from 10mm-30mm.

Indication for use

- It is used in the fixation of Small Fragment Wise-Lock Plates.

