


LCP PROXIMAL FEMUR PLATE (LCP TROCHANTRIC PLATE)



CAT. NO. SS	CAT. NO. TT	LCP PROXIMAL FEMUR PLATE (LCP TROCHANTRIC PLATE)
323.005L	623.005L	LCP TROCHANTRIC PLATE - LEFT, 5H
323.005R	623.005R	LCP TROCHANTRIC PLATE - RIGHT, 5H
323.007L	623.007L	LCP TROCHANTRIC PLATE - LEFT, 7H
323.007R	623.007R	LCP TROCHANTRIC PLATE - RIGHT, 7H
323.009L	623.009L	LCP TROCHANTRIC PLATE - LEFT, 9H
323.009R	623.009R	LCP TROCHANTRIC PLATE - RIGHT, 9H
323.011L	623.011L	LCP TROCHANTRIC PLATE - LEFT, 11H
323.011R	623.011R	LCP TROCHANTRIC PLATE - RIGHT, 11H

Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts LC-LCP Cancellous Screw 6.5mm ,4.5/5.0mm cortical and locking screws.
- Plate available with 5,7,9,11 holes.

Indication for use

- Fractures of the trochanteric region, trochanteric simple, cervicotrochanteric,
- trochanterodiaphyseal, multifragmentary pertrochanteric, intertrochanteric, trochanteric
- reversed or transverse or with additional fracture of the medial cortex
- Fractures of the proximal end of the femur combined with
- ipsilateral shaft fractures
- Metastatic fracture of the proximal femur
- Osteotomies of the proximal femur
- Also for use in fixation of osteopenic