



CAT. NO. SS	CAT. NO. TT	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM
225.016	625.016	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 16 MM
225.018	625.018	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 18 MM
225.020	625.020	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 20 MM
225.022	625.022	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 22 MM
225.024	625.024	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 24 MM
225.026	625.026	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 26 MM
225.028	625.028	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 28 MM
225.030	625.030	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 30 MM
225.032	625.032	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 32 MM
225.034	625.034	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 34 MM
225.036	625.036	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 36 MM
225.038	625.038	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 38 MM
225.040	625.040	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 40 MM
225.042	625.042	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 42 MM
225.044	625.044	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 44 MM
225.046	625.046	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 46 MM
225.048	625.048	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 48 MM
225.050	625.050	LOCKING BOLT (SELF TAPPING) Ø 3.4 MM X 50 MM

Specification

- 3.4mm Locking Bolt is self tapping. Self tapping bolt cuts its own thread while being driven into the bone. It makes a small hole while entering the bone which creates a tight friction fit between the threads. This helps fight vibration loosening and allows the parts to be taken apart if needed.
- Available in both Stainless still & Titanium.
- Available in 16mm to 50mm length with 3.4mm diameter.

Indication for use

- 3.4mm Locking Bolt is used for the fixation of proximal and distal holes of PFNA2, Proximal Femoral Nail (PFN), Antigrade Femoral Nail (AFN), Expert Tibia Nail, Tibia Nail and Femoral Nail.