



CAT. NO. SS	CAT. NO. TT	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM
223.012	623.012	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 12 MM
223.014	623.014	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 14 MM
223.016	623.016	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 16 MM
223.018	623.018	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 18 MM
223.020	623.020	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 20 MM
223.024	623.024	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 24 MM
223.026	623.026	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 26 MM
223.028	623.028	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 28 MM
223.030	623.030	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 30 MM
223.032	623.032	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 32 MM
223.034	623.034	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 34 MM
223.035	623.035	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 35 MM
223.036	623.036	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 36 MM
223.038	623.038	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 38 MM
223.040	623.040	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 40 MM
223.042	623.042	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 42 MM
223.044	623.044	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 44 MM
223.045	623.045	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 45 MM
223.046	623.046	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 46 MM
223.048	623.048	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 48 MM
223.050	623.050	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 50 MM
223.052	623.052	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 52 MM
223.054	623.054	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 54 MM
223.056	623.056	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 56 MM
223.058	623.058	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 58 MM
223.060	623.060	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 60 MM
223.065	623.065	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 65 MM
223.070	623.070	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 70 MM
223.075	623.075	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 75 MM
223.080	623.080	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 80 MM
223.085	623.085	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 85 MM
223.086	623.086	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 86 MM
223.088	623.088	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 88 MM
223.090	623.090	LOCKING BOLT (SELF TAPPING) Ø 4.9 MM X 90 MM

Specification

- 4.9 mm Locking Bolt is self tapping. Self tapping bolt cuts its own thread while being driven into the bone. It makes a small hole while entering the bone which creates a tight friction fit between the threads. This helps fight vibration loosening and allows the parts to be taken apart if needed.
- Available in both Stainless still & Titanium.
- Available in 12mm to 90mm length with 4.9mm diameter.

Indication for use

- 4.9mm Locking Bolt is used for the fixation of proximal and distal holes of PFNA2, Proximal Femoral Nail (PFN), Antigrade Femoral Nail (AFN), Expert Tibia Nail, Tibia Nail and Femoral Nail.