



CAT.NO. SS	CAT. NO. TT	MALLEOLAR SCREWS Ø 4.5MM
216.025	631.025	MALLEOLAR SCREWS Ø 4.5MM X 025 MM
216.030	631.030	MALLEOLAR SCREWS Ø 4.5MM X 030 MM
216.035	631.035	MALLEOLAR SCREWS Ø 4.5MM X 035 MM
216.040	631.040	MALLEOLAR SCREWS Ø 4.5MM X 040 MM
216.045	631.045	MALLEOLAR SCREWS Ø 4.5MM X 045 MM
216.050	631.050	MALLEOLAR SCREWS Ø 4.5MM X 050 MM
216.055	631.055	MALLEOLAR SCREWS Ø 4.5MM X 055 MM
216.060	631.060	MALLEOLAR SCREWS Ø 4.5MM X 060 MM
216.065	631.065	MALLEOLAR SCREWS Ø 4.5MM X 065 MM
216.070	631.070	MALLEOLAR SCREWS Ø 4.5MM X 070 MM
216.075	631.075	MALLEOLAR SCREWS Ø 4.5MM X 075 MM
216.080	631.080	MALLEOLAR SCREWS Ø 4.5MM X 080 MM
216.085	631.085	MALLEOLAR SCREWS Ø 4.5MM X 085 MM
216.090	631.090	MALLEOLAR SCREWS Ø 4.5MM X 090 MM
216.095	631.095	MALLEOLAR SCREWS Ø 4.5MM X 095 MM
216.100	631.100	MALLEOLAR SCREWS Ø 4.5MM X 100 MM

Specification

- Available in 25mm to 100mm length with 4.5mm diameter
- Partially threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- 4.5mm Malleolar screws are indicated for the fixation of large fragments in the cancellous bones

