



CAT. NO.	METAL	Malleolar Screws Ø 3.5mm
216.324	SS	Malleolar Screws Ø 3.5mm x 24 mm
216.326	SS	Malleolar Screws Ø 3.5mm x 26 mm
216.328	SS	Malleolar Screws Ø 3.5mm x 28 mm
216.330	SS	Malleolar Screws Ø 3.5mm x 30 mm
216.332	SS	Malleolar Screws Ø 3.5mm x 32 mm
216.334	SS	Malleolar Screws Ø 3.5mm x 34 mm
216.335	SS	Malleolar Screws Ø 3.5mm x 35 mm
216.336	SS	Malleolar Screws Ø 3.5mm x 36 mm
216.338	SS	Malleolar Screws Ø 3.5mm x 38 mm
216.340	SS	Malleolar Screws Ø 3.5mm x 40 mm
216.342	SS	Malleolar Screws Ø 3.5mm x 42 mm
216.344	SS	Malleolar Screws Ø 3.5mm x 44 mm
216.346	SS	Malleolar Screws Ø 3.5mm x 46 mm
216.348	SS	Malleolar Screws Ø 3.5mm x 48 mm
216.350	SS	Malleolar Screws Ø 3.5mm x 50 mm
216.352	SS	Malleolar Screws Ø 3.5mm x 52mm

Specification

- Available in 20mm to 52mm length with 3.5mm diameter.
- Available in Stainless Steel.

Indication for use

- 4.5mm Malleolar Screw is designed to be used for fracture fixation of small fragments in cancellous bone.
- It is primarily used in inter-fragmental compression bone plates and used in hard cortical where buttress thread form provides enhanced pull-out resistance.