



CAT. NO. SS	CAT. NO. TT	NARROW DCP FOR 4.5 MM SCREWS
276.304	676.304	NARROW DCP FOR 4.5 MM SCREWS - 4 HOLES, L - 87 MM
276.305	676.305	NARROW DCP FOR 4.5 MM SCREWS - 5 HOLES, L - 87 MM
276.306	676.306	NARROW DCP FOR 4.5 MM SCREWS - 6 HOLES, L - 103 MM
276.307	676.307	NARROW DCP FOR 4.5 MM SCREWS - 7 HOLES, L - 119 MM
276.308	676.308	NARROW DCP FOR 4.5 MM SCREWS - 8 HOLES, L - 135 MM
276.309	676.309	NARROW DCP FOR 4.5 MM SCREWS - 9 HOLES, L - 151 MM
276.310	676.310	NARROW DCP FOR 4.5 MM SCREWS - 10 HOLES, L - 167 MM
276.311	676.311	NARROW DCP FOR 4.5 MM SCREWS - 11HOLES, L - 183 MM
276.312	676.312	NARROW DCP FOR 4.5 MM SCREWS - 12 HOLES, L - 199 MM
276.313	676.313	NARROW DCP FOR 4.5 MM SCREWS - 13 HOLES, L - 215 MM
276.314	676.314	NARROW DCP FOR 4.5 MM SCREWS - 14 HOLES, L - 231 MM
276.315	676.315	NARROW DCP FOR 4.5 MM SCREWS - 15 HOLES, L - 247 MM
276.316	676.316	NARROW DCP FOR 4.5 MM SCREWS - 16 HOLES, L - 263 MM

Specification

- Available in both Stainless Steel and Titanium material.
- Plate holes accept 4.5/6.5 mm cortical screws and cancellous screws
- Plates available with 4,5,6,7,8,9,10,11,12,13,14,15 & 16 holes.

Indication for use

4.5mm Narrow Dynamic Compression Plate indicated for fixation of -

- Various long bones, such as the humerus, tibia and femur
- Peri-prosthetic fractures and fractures in osteopenic bone
- Nonunions or malunions in adult patients

