

PELVIC RECONSTRUCTION PLATE

LOW PROFILE (S.S.) A.O TYPE - CURVED 108°



CAT.NO. LVM	CAT.NO. TIT	Pelvic Reconstruction Plate Low Profile (S.S.) A.O Type – Curved 108°
306.106	700.506	PELVIC RECON PLATE 108° Ø 3.5MM - 6 HOLE, CURVED
306.107	700.507	PELVIC RECON PLATE 108° Ø 3.5MM - 7 HOLE, CURVED
306.108	700.508	PELVIC RECON PLATE 108° Ø 3.5MM - 8 HOLE, CURVED
306.110	700.510	PELVIC RECON PLATE 108° Ø 3.5MM - 10 HOLE, CURVED
306.112	700.512	PELVIC RECON PLATE 108° Ø 3.5MM - 12 HOLE, CURVED
306.114	700.514	PELVIC RECON PLATE 108° Ø 3.5MM - 14 HOLE, CURVED
306.116	700.516	PELVIC RECON PLATE 108° Ø 3.5MM - 16 HOLE, CURVED

Specification

- Available in both Titanium & LVM.
- Plate features compression holes that accepts 3.5mm cortical Screws & 4.0mm cancellous screws.
- Plate available with 6,7,8,10,12,14 & 16 holes.

Indication for use

- It is used for pelvic and acetabular reconstruction surgery and fracture fixation of the distal humerus, clavicle and scapula.

