



CAT. NO. LVM	CAT.NO. TT	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM
305.002		PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 2 HOLE, STRAIGHT
305.003		PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 3 HOLE, STRAIGHT
305.004	700.904	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 4 HOLE, STRAIGHT
305.005	700.905	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 5 HOLE, STRAIGHT
305.006	700.906	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 6 HOLE, STRAIGHT
305.007	700.907	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 7 HOLE, STRAIGHT
305.008	700.908	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 8 HOLE, STRAIGHT
305.009	700.909	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 9 HOLE, STRAIGHT
305.010	700.910	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 10 HOLE, STRAIGHT
305.011	700.911	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 11 HOLE, STRAIGHT
305.012	700.912	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 12 HOLE, STRAIGHT
305.013	700.913	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 13 HOLE, STRAIGHT
305.014	700.914	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 14 HOLE, STRAIGHT
305.015	700.915	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 15 HOLE, STRAIGHT
305.016	700.916	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 16 HOLE, STRAIGHT
305.017	700.917	PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 17 HOLE, STRAIGHT
305.018		PELVIC RECON PLATE LOW PROFILE Ø 3.5MM - 18 HOLE, STRAIGHT

Specification

- Available in both Titanium & LVM.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical Screws & 4.0mm cancellous screws.
- Plate available with 2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17 & 18 holes.

Indication for use

- It is used to fix Pelvic, clavicle and calcaneal fractures.