

PELVIC RECONSTRUCTION PLATE

LOW PROFILE (S.S.) A.O TYPE - CURVED 88°



CAT. NO. LVM	CAT. NO. TIT	Pelvic Reconstruction Plate Low Profile (S.S.) A.O Type - Curved 88°
305.106	700.486	PELVIC RECON PLATE 88° Ø 3.5MM - 6 HOLE, CURVED
305.107		PELVIC RECON PLATE 88° Ø 3.5MM - 7 HOLE, CURVED
305.108	700.488	PELVIC RECON PLATE 88° Ø 3.5MM - 8 HOLE, CURVED
305.110	700.490	PELVIC RECON PLATE 88° Ø 3.5MM - 10 HOLE, CURVED
305.112	700.492	PELVIC RECON PLATE 88° Ø 3.5MM - 12 HOLE, CURVED
305.114	700.494	PELVIC RECON PLATE 88° Ø 3.5MM - 14 HOLE, CURVED
305.115		PELVIC RECON PLATE 88° Ø 3.5MM - 15 HOLE, CURVED
305.116	700.496	PELVIC RECON PLATE 88° Ø 3.5MM - 16 HOLE, CURVED

Specification

- Available in LVM.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plate available with 2,3,4,5,6,7,8,9,10,12,14,16 & 18 holes.

Indication for use

- It is used to fix Pelvic, clavicle and calcaneal fractures.

