

# **PELVIC RECONSTRUCTION** **PLATE LOW PROFILE (S.S.) A.O TYPE - J PLATE**



CAT. NO. LVM	CAT. NO. TT	Pelvic Reconstruction Plate Low Profile (S.S.) A.O Type - J Plate
305.210L	700.823L	PELVIC J PLATE Ø 3.5MM - 10 HOLE, LEFT
305.210R	700.823R	PELVIC J PLATE Ø 3.5MM - 10 HOLE, RIGHT
305.212L	700.825L	PELVIC J PLATE Ø 3.5MM - 12 HOLE, LEFT
305.212R	700.825R	PELVIC J PLATE Ø 3.5MM - 12 HOLE, RIGHT
305.214L	700.827L	PELVIC J PLATE Ø 3.5MM - 14 HOLE, LEFT
305.214R	700.827R	PELVIC J PLATE Ø 3.5MM - 14 HOLE, RIGHT
305.216L	700.829L	PELVIC J PLATE Ø 3.5MM - 16 HOLE, LEFT
305.216R	700.829R	PELVIC J PLATE Ø 3.5MM - 16 HOLE, RIGHT

## **Specification**

- Available in both Titanium & LVM.
- Plate features compression hole that accepts 3.5mm cortical Screws & 4.0mm cancellous screws.
- Plate available with 10,12,14 & 16 holes.

## **Indication for use**

- Pelvic Reconstruction Plate System are intended for pelvic and acetabular reconstructive surgery.
- Additional implant tray for Low Profile 3.5 Reconstruction J-Plate Set.
- J-Plates are pre-contoured to fit the pelvic brim.

