



CAT. NO.	METAL	Posts Female
382.001	SS	Posts Female - 01 Holes
382.002	SS	Posts Female - 02 Holes
382.003	SS	Posts Female - 03 Holes
382.004	SS	Posts Female - 04 Holes

Specification

- Posts Female available in Stainless steel.
- It is available in 1,2,3 & 4 holes.

Indication for use

- Posts Female is indicated as an external fixator, it gives relative stability. as pins are inserted across different planes in a multiplanar fixation, the construct provides great stability.