



CAT. NO.	METAL	PostsMale
383.001	SS	Posts Male - 01 Holes
383.002	SS	Posts Male - 02 Holes
383.003	SS	Posts Male - 03 Holes
383.004	SS	Posts Male - 04 Holes

## Specification

- Posts Male available in Stainless steel.
- It is available in 1,2,3 & 4 holes.

## Indication for use

- Posts Male is indicated as an external fixator, it gives relative stability. as pins are inserted across different planes in a multiplanar fixation, the construct provides great stability.