



| CAT.NO. SS | CAT. NO. TT | SMALL DCP FOR 3.5 MM SCREWS |
|---------------|----------------|--|
| 276.205 | 676.20 5 | SMALL DCP FOR 3.5 MM SCREWS - 05 HOLES, L - 62 MM |
| 276.206 | 676.20 6 | SMALL DCP FOR 3.5 MM SCREWS - 06 HOLES, L - 77 MM |
| 276.207 | 676.20 7 | SMALL DCP FOR 3.5 MM SCREWS - 07 HOLES, L - 86 MM |
| 276.208 | 676.20 8 | SMALL DCP FOR 3.5 MM SCREWS - 08 HOLES, L - 98 MM |
| 276.209 | 676.20 9 | SMALL DCP FOR 3.5 MM SCREWS - 09 HOLES, L - 110 MM |
| 276.210 | 676.21 0 | SMALL DCP FOR 3.5 MM SCREWS - 10 HOLES, L - 122 MM |
| 276.211 | 676.21 1 | SMALL DCP FOR 3.5 MM SCREWS - 11 HOLES, |
| 276.212 | 676.21 2 | SMALL DCP FOR 3.5 MM SCREWS - 12 HOLES, |
| 276.213 | 676.21 3 | SMALL DCP FOR 3.5 MM SCREWS - 13 HOLES, |
| 276.214 | 676.21 4 | SMALL DCP FOR 3.5 MM SCREWS - 14 HOLES, |

Specification

- Available in both Stainless Steel and Titanium material.
- The oval shape holes are designed to create dynamic pressure between the fractured bone fragments.
- All holes of the plate support 3.5mm cortical & 4.0mm cancellous screws.

Indication for use

- It is used to fix fractures and non-unions of the fibula, Distal tibia, Pelvis, Ulna, Radius, Humorous, Olecranon, Scapula & Clavicle.

