

## SMALL T- PLATE- OBLIQUE ANGLE (LEFT & RIGHT COMBINED)



CAT. NO. SS	Small T - Plate- Oblique Angle (Left & Right Combined)
295.030	Small T - Plate- Oblique Angle - 03 Hole, L - 52mm
295.040	Small T - Plate- Oblique Angle - 04 Hole, L - 63mm
295.050	Small T - Plate- Oblique Angle - 05 Hole, L - 74mm
295.060	Small T - Plate- Oblique Angle - 06 Hole,
295.070	Small T - Plate- Oblique Angle - 07 Hole,
295.080	Small T - Plate- Oblique Angle - 08 Hole,

### Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5,6 & 7 Holes.
- Plate allows 3.5 mm cortical and 4.0 mm cancellous Screws.

### Indication for use

- These plates are indicated for Extra-articular, Partial articular and Complete articular fractures of Distal Radius.
- These plates are designed separately designed for left and right hand. All holes of the plate allow 3.5 mm cortical and 4.0 mm cancellous Screws.

