



CAT. NO. SS	CAT. NO. TT	SMALL T - PLATE - RIGHT ANGLE 3.5mm Screw
295.033	699.033	SMALL T - PLATE - RIGHT ANGLE - 03/3 (SHAFT/ HEAD) HOLE, L - 50MM
295.034	699.034	SMALL T - PLATE - RIGHT ANGLE - 03/4 (SHAFT/ HEAD) HOLE, L - 50MM
295.043	699.043	SMALL T - PLATE - RIGHT ANGLE - 04/3 (SHAFT/ HEAD) HOLE, L - 56MM
295.044	699.044	SMALL T - PLATE - RIGHT ANGLE - 04/4 (SHAFT/ HEAD) HOLE, L - 56MM
295.053	699.053	SMALL T - PLATE - RIGHT ANGLE - 05/3 (SHAFT/ HEAD) HOLE, L - 67MM
295.054	699.054	SMALL T - PLATE - RIGHT ANGLE - 05/4 (SHAFT/ HEAD) HOLE, L - 67MM
295.063	699.063	SMALL T - PLATE - RIGHT ANGLE - 06/3 (SHAFT/ HEAD) HOLE, L - 78MM
295.064	699.064	SMALL T - PLATE - RIGHT ANGLE - 06/4 (SHAFT/ HEAD) HOLE, L - 78MM
295.074	699.074	SMALL T - PLATE - RIGHT ANGLE - 07/4 (SHAFT/ HEAD) HOLE,
295.084	699.084	SMALL T - PLATE - RIGHT ANGLE - 08/4 (SHAFT/ HEAD) HOLE,

Specification

- Available in both Titanium & Stainless Steel.
- All holes of the plate support 3.5mm cortical & 4.0mm cancellous screws.

Indication for use

- These plates are indicated for the extra-articular, partial articular and complete articular fractures of distal radius(left & right radius).