



CAT. NO. LVM	CAT. NO. TT	SPRING PLATE 3.5
305.251	700.836	SPRING PLATE 3.5 -1H LENGTH 19.5MM
305.252	700.837	SPRING PLATE 3.5 -2H LENGTH 31.5MM
305.253	700.838	SPRING PLATE 3.5 -3H LENGTH 43.5MM

Specification

- Available in both in LVM & Titanium
- This plates are part of the 3.5 mm Low Profile Pelvic System, which offers implants and instruments for pelvic and acetabular reconstructive surgery.
- Available in 1,2 & 3 Holes .

Indication for use

- Spring Plates used for Reduction and fixation of small bone fragments.
- Spring Plate is intended for pelvic and acetabular reconstructive surgery and fracture fixation of the distal fibula.

