



CAT. NO. SS	CAT. NO. TT	ST. RECON PLATE - 1.6 MM X 3.5MM
292.005	709.005	ST. RECON PLATE - 1.6 MM X 3.5MM - 05 HOLES
292.006	709.006	ST. RECON PLATE - 1.6 MM X 3.5MM - 06 HOLES
292.007	709.007	ST. RECON PLATE - 1.6 MM X 3.5MM - 07 HOLES
292.008	709.008	ST. RECON PLATE - 1.6 MM X 3.5MM - 08 HOLES
292.009	709.009	ST. RECON PLATE - 1.6 MM X 3.5MM - 09 HOLES
292.010	709.010	ST. RECON PLATE - 1.6 MM X 3.5MM - 10 HOLES
292.011	709.011	ST. RECON PLATE - 1.6 MM X 3.5MM - 11 HOLES
292.012	709.012	ST. RECON PLATE - 1.6 MM X 3.5MM - 12 HOLES
292.013	709.013	ST. RECON PLATE - 1.6 MM X 3.5MM - 13 HOLES
292.014	709.014	ST. RECON PLATE - 1.6 MM X 3.5MM - 14 HOLES
292.015	709.015	ST. RECON PLATE - 1.6 MM X 3.5MM - 15 HOLES
292.016	709.016	ST. RECON PLATE - 1.6 MM X 3.5MM - 16 HOLES
292.018	709.018	ST. RECON PLATE - 1.6 MM X 3.5MM - 18 HOLES

Specification

- These plates are used with 3.5mm cortical & 4.0mm cancellous screws.
- Plates bears side cuts which makes them more flexible.
- Available in both Titanium & Stainless Steel.

Indication for use

- Reconstruction Plate Straight is used to fix Diaphysis Noncomminuted (Spiral, Oblique or Transverse) and Diaphysis Wedge (Spiral, Bending or Comminuted) fractures of metacarpals and metatarsals.

