




CAT. NO. SS	CAT. NO. TT	ST. RECON PLATE FOR 3.5 MM SCREWS
292.303	709.303	ST. RECON PLATE FOR 3.5 MM SCREWS - 03 HOLES, L - 34MM
292.304	709.304	ST. RECON PLATE FOR 3.5 MM SCREWS - 04 HOLES, L - 46MM
292.305	709.305	ST. RECON PLATE FOR 3.5 MM SCREWS - 05 HOLES, L - 58MM
292.306	709.306	ST. RECON PLATE FOR 3.5 MM SCREWS - 06 HOLES, L - 70MM
292.307	709.307	ST. RECON PLATE FOR 3.5 MM SCREWS - 07 HOLES, L - 82MM
292.308	709.308	ST. RECON PLATE FOR 3.5 MM SCREWS - 08 HOLES, L - 94MM
292.309	709.309	ST. RECON PLATE FOR 3.5 MM SCREWS - 09 HOLES, L - 106MM
292.310	709.310	ST. RECON PLATE FOR 3.5 MM SCREWS - 10 HOLES, L - 118MM
292.311	709.311	ST. RECON PLATE FOR 3.5 MM SCREWS - 11 HOLES, L - 130MM
292.312	709.312	ST. RECON PLATE FOR 3.5 MM SCREWS - 12 HOLES, L - 142MM
292.313	709.313	ST. RECON PLATE FOR 3.5 MM SCREWS - 13 HOLES, L - 154MM
292.314	709.314	ST. RECON PLATE FOR 3.5 MM SCREWS - 14 HOLES, L - 166MM
292.315	709.315	ST. RECON PLATE FOR 3.5 MM SCREWS - 15 HOLES, L - 178MM
292.316	709.316	ST. RECON PLATE FOR 3.5 MM SCREWS - 16 HOLES, L - 190MM
292.317	709.317	ST. RECON PLATE FOR 3.5 MM SCREWS - 17 HOLES, L - 202MM
292.318	709.318	ST. RECON PLATE FOR 3.5 MM SCREWS - 18 HOLES, L - 214MM

Specification

- These plates are used with 3.5mm cortical & 4.0mm cancellous screws.
- Plates bears side cuts which makes them more flexible.
- Available in both Titanium & Stainless Steel.

Indication for use

- 3.5mm Reconstruction Plate, Straight is used fix clavicle and calcaneal fractures.

