



CAT. NO.	METAL	Schanz Screws- (Front Threaded)
232.250	SS	Schanz Screws Ø 2.5 mm x 14th x 100mm
232.251	SS	Schanz Screws Ø 3.5 mm x 16th x 125mm
232.252	SS	Schanz Screws Ø 3.5 mm x 16th x 150mm
232.265	SS	Schanz Screws Ø 3.5 mm x 32th x 150mm
232.280	SS	Schanz Screws Ø 3.5 mm x 32th x 225mm
232.406	SS	Schanz Screws Ø 4.5 x16th x 150mm
232.407	SS	Schanz Screws Ø 4.5 x16th x 175mm
232.408	SS	Schanz Screws Ø 4.5 x16th x 200mm
232.409	SS	Schanz Screws Ø 4.5 x16th x 225mm
232.412	SS	SchanzScrews Ø 4.5 x18th x 150mm
232.413	SS	Schanz Screws Ø 4.5 x18th x 175mm
232.414	SS	Schanz Screws Ø 4.5 x18th x 200mm
232.415	SS	Schanz Screws Ø 4.5 x18th x 225mm
232.417	SS	Schanz Screws Ø 4.5 mm x 25th x 150mm
232.418	SS	Schanz Screws Ø 4.5 mm x 25th x 175mm
232.419	SS	Schanz Screws Ø 4.5 mm x 25th x 200mm
232.420	SS	Schanz Screws Ø 4.5 mm x 25th x 250mm
232.421	SS	Schanz Screws Ø 4.5 mm x 32th x 150mm
232.422	SS	Schanz Screws Ø 4.5 mm x 32th x 175mm
232.423	SS	Schanz Screws Ø 4.5 mm x 32th x 200mm
232.424	SS	Schanz Screws Ø 4.5 mm x 32th x 225mm
232.425	SS	Schanz Screws Ø 4.5 mm x 40th x 175mm
232.426	SS	Schanz Screws Ø 4.5 mm x 40th x 200mm
232.427	SS	Schanz Screws Ø 4.5 mm x 40th x 225mm
232.428	SS	Schanz Screws Ø 4.5 mm x 40th x 150mm
232.455	SS	Schanz Screws Ø 4.5 mm x 25th x 225mm
232.503	SS	Schanz Screws Ø 5.0 mm x 18th x 200mm
232.511	SS	Schanz Screws Ø 5.0 mm x 25th x 150mm
232.523	SS	Schanz Screws Ø 5.0 mm x 32th x 200mm
232.524	SS	Schanz Screws Ø 5.0 mm x 32th x 225mm
232.528	SS	Schanz Screws Ø 5.0 mm x 40th x 200mm
232.535	SS	Schanz Screws Ø 5.0 mm x 50th x 200mm
232.536	SS	Schanz Screws Ø 5.0 mm x 50th x 225mm
232.603	SS	Schanz Screws Ø 6.0 mm x 32th x 200mm
232.604	SS	Schanz Screws Ø 6.0 mm x 32th x 225mm
232.614	SS	Schanz Screws Ø 6.0 mm x 50th x 200mm
232.615	SS	Schanz Screws Ø 6.0 mm x 50th x 225mm
232.631	SS	Schanz Screws Ø 6.5mm x 40th x 225mm
323.712	SS	Schanz Screws Ø 4.0 x25th x 175mm
323.714	SS	Schanz Screws Ø 4.0 x40th x 175mm
357.326	SS	Schanz Screws Ø 6.0 mm x 32th x 130mm

Specification

- Schanz Screws - (Front Threaded) available in Stainless steel.
- It is available in 100,125,130,150,175,200,225 & 250mm length with 2.5,3.5,4.5,5.0 & 6.0mm Dia.

Indication for use

- Schanz Screw is intended for use with an external fixation system for fracture fixation (open or closed); pseudoarthrosis or nonunion of long bones; limb lengthening by epiphyseal or metaphyseal distraction; correction of bony or soft tissue deformity; correction of segmental bony or soft tissue deformity; correction of segmental bony or soft tissue defects; and joint arthrodesis.

