

# SCHANZ SCREWS SELF DRILLING AND SELF TAPPING



CAT. NO.	METAL	Schanz Screws Self drilling and Self tapping
232.617	SS	Schanz Screws Self drilling and Self tapping Ø 5.0 mm x 75th x 200mm
232.618	SS	Schanz Screws Self drilling and Self tapping Ø 4.0 mm x 40th x 150mm

## Specification

- Schanz Screws Self drilling and Self tapping available in Stainless steel.
- It is available in 150 & 200mm length with 4.0 & 5.0mm Dia.

## Indication for use

- Schanz Screws self drilling and self tapping has been specifically designed to optimise the bone/pin interface to reduce the occurrence of pin-tract related complications in external fixation.

