



| CAT. NO. | METAL | Steinmann Pin 3.0 mm |
|----------|-------|------------------------------|
| 235.150 | SS | Steinmann Pin 3.0 mm x 150mm |
| 235.175 | SS | Steinmann Pin 3.0 mm x 175mm |
| 235.200 | SS | Steinmann Pin 3.0 mm x 200mm |
| 235.225 | SS | Steinmann Pin 3.0 mm x 225mm |
| 235.230 | SS | Steinmann Pin 3.0 mm x 300mm |

Specification

- Steinmann Pin available in Stainless steel.
- Available length 150,175,200,225 & 300mm with 3.0mm dia.

Indication for use

- Steinmann Pin is indicated for use in the fixation of bone fractures, bone reconstruction, and as guide pins for insertion of other implant.

