



CAT. NO.	METAL	Steinmann Pin 35 mm
236.125	SS	Steinmann Pin 3.5 mm x 125mm
236.150	SS	Steinmann Pin 3.5 mm x 150mm
236.175	SS	Steinmann Pin 3.5 mm x 175mm
236.200	SS	Steinmann Pin 3.5 mm x 200mm
236.225	SS	Steinmann Pin 3.5 mm x 225mm
236.250	SS	Steinmann Pin 3.5 mm x 250mm
236.260	SS	Steinmann Pin 3.5 mm x 300mm

Specification

- Steinmann Pin 3.5MM available in Stainless steel.
- Available length 125,150,175,200,225,250 & 300mm with 3.5mm dia.

Indication for use

- Steinmann Pin is indicated for use in the fixation of bone fractures, bone reconstruction, and as guide pins for insertion of other implant.

