



CAT. NO.	METAL	Steinmann Pin 4.0 mm
237.145	SS	Steinmann Pin 4.0 mm x 125mm
237.150	SS	Steinmann Pin 4.0 mm x 150mm
237.175	SS	Steinmann Pin 4.0 mm x 175mm
237.180	SS	Steinmann Pin 4.0 mm x 180mm
237.200	SS	Steinmann Pin 4.0 mm x 200mm
237.225	SS	Steinmann Pin 4.0 mm x 225mm
237.250	SS	Steinmann Pin 4.0 mm x 250mm
237.260	SS	Steinmann Pin 4.0 mm x 300mm

Specification

- Steinmann Pin 4.0mm available in Stainless steel.
- Available length 125,150,175,180,200,225,250 & 300mm with 4.0mm dia.

Indication for use

- Steinmann Pin is indicated for use in the fixation of bone fractures, bone reconstruction, and as guide pins for insertion of other implant.

