



CAT. NO.	METAL	Steinmann Pin 4.5 mm
238.150	SS	Steinmann Pin Ø 4.5 mm x 150mm
238.175	SS	Steinmann Pin Ø 4.5 mm x 175mm
238.200	SS	Steinmann Pin Ø 4.5 mm x 200mm
238.225	SS	Steinmann Pin Ø 4.5 mm x 225mm
238.250	SS	Steinmann Pin Ø 4.5 mm x 250mm
238.260	SS	Steinmann Pin Ø 4.5 mm x 300mm

Specification

- Steinmann Pin 4.5mm available in Stainless steel.
- Available length 150,175,200,225,250 & 300mm with 4.5mm dia.

Indication for use

- Steinmann Pin is indicated for use in the fixation of bone fractures, bone reconstruction, and as guide pins for insertion of other implant.

