



CAT. NO.	METAL	Steinmann Pin 5.0 mm
239.150	SS	Steinmann Pin Ø 5.0 mm x 150mm
239.175	SS	Steinmann Pin Ø 5.0 mm x 175mm
239.200	SS	Steinmann Pin Ø 5.0 mm x 200mm
239.225	SS	Steinmann Pin Ø 5.0 mm x 225mm

## Specification

- Steinmann Pin 5.0mm available in Stainless steel.
- Available length 150,175,200,225,250 & 300mm with 4.5mm dia.

## Indication for use

- Steinmann Pin is indicated for use in the fixation of bone fractures, bone reconstruction, and as guide pins for insertion of other implant.

