



CAT. NO.SS	Steinmann Pin- Trocar Tip (Sterilised)
238.225S	Steinmann Pin Ø 4.5 mm x 225mm(Sterilised)
239.225S	Steinmann Pin Ø 5.0 mm x 225mm(Sterilised)

Specification

- Steinmann Pin - Trocar Tip (Sterilised) available in Stainless steel.
- Available length 225mm with 4.5 & 5.0mm dia.

Indication for use

- Steinman pin used for skeletal traction and fixation of bone fragments.

