


STRAIGHT (PELVIC) RECONSTRUCTION PLATES- 4.5MM SCREWS



CAT. NO. SS	Straight (Pelvic) Reconstruction Plates - 4.5mm Screws
292.408	St. Recon Plate for 4.5 mm Screws - 08 Holes, L - 94mm
292.409	St. Recon Plate for 4.5 mm Screws - 09 Holes, L - 106mm
292.410	St. Recon Plate for 4.5 mm Screws - 10 Holes, L - 118mm
292.411	St. Recon Plate for 4.5 mm Screws - 11 Holes, L - 130mm
292.412	St. Recon Plate for 4.5 mm Screws - 12 Holes, L - 142mm
292.413	St. Recon Plate for 4.5 mm Screws - 13 Holes, L - 154mm
292.414	St. Recon Plate for 4.5 mm Screws - 14 Holes, L - 166mm

Specification

- Plates available in LVM.
- Plates available with 8,9,10,11,12,13 & 14 holes.
- Plate accepts 4.5mm Cortical screw.

Indication for use

- These plates are primarily indicated for the reconstruction of the pelvis.
- It is also indicated for condylar, metaphyseal, or diaphyseal fractures of small and long bones.

