



| CAT. NO. SS | T- Buttress Plates for 4.5mm Screws |
|-------------|--|
| 296.003 | T- Buttress Plates - 03 Holes, L - 66 mm |
| 296.00 4 | T- Buttress Plates - 04 Holes, L - 81 mm |
| 296.00 5 | T- Buttress Plates - 05 Holes, L - 96 mm |
| 296.00 6 | T- Buttress Plates - 06 Holes, L - 112mm |
| 296.00 7 | T- Buttress Plates - 07 Holes, L - 128mm |
| 296.00 8 | T- Buttress Plates - 08 Holes, L - 144mm |
| 296.00 9 | T- Buttress Plates - 09 Holes, L - 144mm |
| 296.010 | T- Buttress Plates - 10 Holes, L - 176mm |

Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5,6,7,8,9 & 10 Holes.
- Plates accepts 4.5mm cortical screws.

Indication for use

- T-Buttress Plate for the lateral plateau fracture of the tibia. It is designed to be used as a buttress plate in areas of thin metaphyseal cortical bone or to prevent the collapsing during the treatment of cancellous defects in that area. Prior to the surgery the plate may be bent to accurately match the anatomy of the segment of bone to be fixed. The long oval hole is provided for temporary fixation of the plate prior to final adjustment and fixation.

