



CAT. NO. SS	T- Plate- 4.5MM
297.003	T- Plate- 03 Holes, L - 68 mm
297.004	T- Plate- 04 Holes, L - 84 mm
297.005	T- Plate- 05 Holes, L - 100mm
297.006	T- Plate- 06 Holes, L - 116mm
297.007	T- Plate- 07 Holes, L - 132mm
297.008	T- Plate- 08 Holes, L - 148mm
297.009	T- Plate- 09 Holes, L - 164mm
297.010	T- Plate- 10 Holes, L - 180mm

Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5,6,7,8,9 & 10 holes.
- Plate accepts 4.5mm Cortical screw.

Indication for use

- It is used in the distal tibia, medial distal femur and proximal humerus fractures.
- It is designed to use in metaphyseal regions of bone.

