



| CAT. NO. | METAL | Tapered Schanz Screw |
|----------|-------|--|
| 357.020 | SS | Tapered Schanz Screw Ø 4.5mm x 20th x 100 mm |
| 357.030 | SS | Tapered Schanz Screw Ø 4.5mm x 30th x 100 mm |
| 357.120 | SS | Tapered Schanz Screw Ø 4.5mm x 40th x 200 mm |
| 357.130 | SS | Tapered Schanz Screw Ø 4.5mm x 30th x 120 mm |
| 357.139 | SS | Tapered Schanz Screw Ø 4.5mm x 40th x 100 mm |
| 357.140 | SS | Tapered Schanz Screw Ø 4.5mm x 40th x 120 mm |
| 357.225 | SS | Tapered Schanz Screw Ø 4.5mm x 40th x 150 mm |
| 357.345 | SS | Tapered Schanz Screw Ø 4.5mm x 40th x 170 mm |

Specification

- Tapered Schanz Screw available in Stainless steel.
- Available length 100,120,150,170 & 200mm with 4.5mm dia.

Indication for use

- Schanz Screw is intended for use with an external fixation system for fracture fixation (open or closed); pseudoarthrosis or nonunion of long bones; limb lengthening by epiphyseal or metaphyseal distraction; correction of bony or soft tissue deformity; correction of segmental bony or soft tissue deformity; correction of segmental bony or soft tissue defects; and joint arthrodesis

