



CAT. NO.	METAL	Threaded Socket
390.020	SS	Threaded Socket - 20 mm
390.030	SS	Threaded Socket - 30 mm
390.040	SS	Threaded Socket - 40 mm
390.050	SS	Threaded Socket - 60 mm
390.100	SS	Square Socket

Specification

- Threaded Socket available in Stainless steel.
- It is available in 20,30,40 & 60mm length.

Indication for use

- Threaded Socket device is a specialized form of external fixator, a circular fixator, modular in construction. Stainless steel / Carbon Fiber rings are fixed to the bone via stainless heavy-gauge wire (called “pins” or Kirschner wires). The rings are connected to each other with threaded rods attached through adjustable nuts.

