



CAT. NO.	METAL	U- Type Staples (Epiphyseal)
346.112	SS	U- Type Staples (Epiphyseal) Ø 1.5mm x 12 mm
346.118	SS	U- Type Staples (Epiphyseal) Ø 1.5mm x 18 mm
346.136	SS	U- Type Staples (Epiphyseal) Ø 1.5mm x 36 mm
346.212	SS	U- Type Staples (Epiphyseal) Ø 2.0mm x 12 mm
346.224	SS	U- Type Staples (Epiphyseal) Ø 2.0mm x 24 mm
346.236	SS	U- Type Staples (Epiphyseal) Ø 2.0mm x 36 mm

Specification

- U - Type Staples (Epiphyseal) available in Stainless steel.
- Available length 12,18,24 & 36mm with 1.5mm and 2.0mm dia.

Indication for use

- Epiphyseal Staple is used for correcting angular deformity at the knee. Two or more Epiphyseal Staple should be used on each side of the epiphysis.

