



CAT. NO. SS	CAT. NO. TT	PFNA II - BLADES (ANTIROTATION SCREW)
268.708	668.708	PFNA II BLADE 60.0 MM
268.709	668.709	PFNA II BLADE 65.0 MM
268.710	668.710	PFNA II BLADE 70.0 MM
268.711	668.711	PFNA II BLADE 75.0 MM
268.712	668.712	PFNA II BLADE 80.0 MM
268.713	668.713	PFNA II BLADE 85.0 MM
268.714	668.714	PFNA II BLADE 90.0 MM
268.715	668.715	PFNA II BLADE 95.0 MM
268.716	668.716	PFNA II BLADE 100.0 MM
268.717	668.717	PFNA II BLADE 105.0 MM
268.718	668.718	PFNA II BLADE 110.0 MM
268.719	668.719	PFNA II BLADE 115.0 MM
268.720	668.720	PFNA II BLADE 120.0 MM

Specification

- Available in both Stainless still & Titanium.
- Available in 60 to 120mm diameter.
- The PFNA-II blade is automatically locked to prevent rotation of the blade and femoral head

Indication for use

- Pertrochanteric fractures
- Intertrochanteric fractures
- High subtrochanteric fractures
- Low and extended subtrochanteric fractures
- Combination fractures (in the proximal femur)
- Pathological fractures

