

CAT. NO. SS	SQUARE NAIL FOR RADIUS
248.015	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 15 CM
248.016	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 16 CM
248.017	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 17 CM
248.018	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 18 CM
248.019	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 19 CM
248.020	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 20 CM
248.021	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 21 CM
248.022	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 22 CM
248.023	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 23 CM
248.024	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 24 CM
248.025	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 25 CM
248.026	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 26 CM
248.027	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 27 CM
248.028	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 28 CM
248.029	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 29 CM
248.030	SQUARE NAIL FOR RADIUS Ø 2.0 MM X 30 CM
248.115	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 15 CM
248.116	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 16 CM
248.117	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 17 CM
248.118	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 18 CM
248.119	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 19 CM
248.120	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 20 CM
248.121	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 21 CM
248.122	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 22 CM
248.123	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 23 CM
248.124	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 24 CM
248.125	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 25 CM
248.126	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 26 CM
248.127	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 27 CM
248.128	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 28 CM
248.129	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 29 CM
248.130	SQUARE NAIL FOR RADIUS Ø 2.5 MM X 30 CM
248.215	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 15 CM
248.216	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 16 CM
248.217	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 17 CM
248.218	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 18 CM
248.219	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 19 CM
248.220	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 20 CM
248.221	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 21 CM
248.222	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 22 CM
248.223	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 23 CM
248.224	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 24 CM
248.225	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 25 CM
248.226	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 26 CM
248.227	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 27 CM
248.228	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 28 CM
248.229	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 29 CM
248.230	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 30 CM
248.236	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 36 CM
248.237	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 37 CM
248.238	SQUARE NAIL FOR RADIUS Ø 3.0 MM X 38 CM
248.315	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 15 CM
248.316	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 16 CM
248.317	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 17 CM
248.318	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 18 CM
248.319	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 19 CM
248.320	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 20 CM
248.321	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 21 CM
248.322	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 22 CM
248.323	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 23 CM
248.324	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 24 CM

248.325	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 25 CM
248.326	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 26 CM
248.327	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 27 CM
248.328	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 28 CM
248.329	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 29 CM
248.330	SQUARE NAIL FOR RADIUS Ø 3.5 MM X 30 CM
248.416	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 16 CM
248.417	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 17 CM
248.418	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 18 CM
248.419	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 19 CM
248.420	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 20 CM
248.421	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 21 CM
248.422	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 22 CM
248.423	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 23 CM
248.424	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 24 CM
248.425	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 25 CM
248.426	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 26 CM
248.427	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 27 CM
248.428	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 28 CM
248.429	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 29 CM
248.430	SQUARE NAIL FOR RADIUS Ø 4.0 MM X 30 CM

Specification

- Available in both Stainless steel
- Available in 15cm to 30cm length with 2.0, 2.5, 3.0, 3.5 & 4.0mm diameter.

Indication for use

- Square Nails for Radius used to treat Diaphyseal fractures of Radius Bone.

The square shaft prevents rotation and stabilizes the fixation.