

| CAT. NO.<br>SS | SQUARE NAIL FOR ULNA                  |
|----------------|---------------------------------------|
| 249.015        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 15 CM |
| 249.016        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 16 CM |
| 249.017        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 17 CM |
| 249.018        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 18 CM |
| 249.019        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 19 CM |
| 249.020        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 20 CM |
| 249.021        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 21 CM |
| 249.022        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 22 CM |
| 249.023        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 23 CM |
| 249.024        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 24 CM |
| 249.025        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 25 CM |
| 249.026        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 26 CM |
| 249.027        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 27 CM |
| 249.028        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 28 CM |
| 249.029        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 29 CM |
| 249.030        | SQUARE NAIL FOR ULNA Ø 2.0 MM X 30 CM |
| 249.115        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 15 CM |
| 249.116        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 16 CM |
| 249.117        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 17 CM |
| 249.118        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 18 CM |
| 249.119        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 19 CM |
| 249.120        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 20 CM |
| 249.121        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 21 CM |
| 249.122        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 22 CM |
| 249.123        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 23 CM |
| 249.124        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 24 CM |
| 249.125        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 25 CM |
| 249.126        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 26 CM |
| 249.127        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 27 CM |
| 249.128        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 28 CM |
| 249.129        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 29 CM |
| 249.130        | SQUARE NAIL FOR ULNA Ø 2.5 MM X 30 CM |
| 249.215        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 15 CM |
| 249.216        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 16 CM |
| 249.217        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 17 CM |
| 249.218        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 18 CM |
| 249.219        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 19 CM |
| 249.220        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 20 CM |
| 249.221        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 21 CM |
| 249.222        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 22 CM |
| 249.223        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 23 CM |
| 249.224        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 24 CM |
| 249.225        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 25 CM |
| 249.226        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 26 CM |
| 249.227        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 27 CM |
| 249.228        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 28 CM |
| 249.229        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 29 CM |
| 249.230        | SQUARE NAIL FOR ULNA Ø 3.0 MM X 30 CM |
| 249.315        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 15 CM |
| 249.316        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 16 CM |
| 249.317        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 17 CM |
| 249.318        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 18 CM |
| 249.319        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 19 CM |
| 249.320        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 20 CM |
| 249.321        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 21 CM |
| 249.322        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 22 CM |
| 249.323        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 23 CM |
| 249.324        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 24 CM |
| 249.325        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 25 CM |
| 249.326        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 26 CM |
| 249.327        | SQUARE NAIL FOR ULNA Ø 3.5 MM X 27 CM |

|         |                                       |
|---------|---------------------------------------|
| 249.328 | SQUARE NAIL FOR ULNA Ø 3.5 MM X 28 CM |
| 249.329 | SQUARE NAIL FOR ULNA Ø 3.5 MM X 29 CM |
| 249.330 | SQUARE NAIL FOR ULNA Ø 3.5 MM X 30 CM |
| 249.416 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 16 CM |
| 249.417 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 17 CM |
| 249.418 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 18 CM |
| 249.419 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 19 CM |
| 249.420 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 20 CM |
| 249.421 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 21 CM |
| 249.422 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 22 CM |
| 249.423 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 23 CM |
| 249.424 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 24 CM |
| 249.425 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 25 CM |
| 249.426 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 26 CM |
| 249.427 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 27 CM |
| 249.428 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 28 CM |
| 249.429 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 29 CM |
| 249.430 | SQUARE NAIL FOR ULNA Ø 4.0 MM X 30 CM |

## Specification

- Available in Stainless steel.
- Available in 15cm to 30cm length with 2.0, 2.5, 3.0, 3.5 & 4.0mm diameter.

## Indication for use

- Square Nails for Ulna used to treat Diaphyseal fractures of Ulna bone.

The square shaft prevents rotation and stabilizes the fixation.