



CAT. NO. SS	CAT. NO. TT	INTER LOCKING ARTHODEHSIS NAIL
256.070	656.070	INTER LOCKING ARTHODEHSIS NAIL 09 MM /10MM X 70 CM
256.075	656.075	INTER LOCKING ARTHODEHSIS NAIL 09 MM /10MM X 75 CM
256.080	656.080	INTER LOCKING ARTHODEHSIS NAIL 09 MM /10MM X 80 CM

## Specification

- Available in both Stainless steel & Titanium.
- Available in 70cm to 80cm length with 9.0 & 10mm diameter.

## Indication for use

- It is used for intramedullary fixation for knee joint fusion.

