

LONG PROXIMAL FEMORAL NAILS 130*/135* (LONG P.F.N.)



CAT. NO. SS	CAT. NO. TT	Long Proximal Femoral Nails 130*/135* (Long P.F.N.)
267.030L	667.030L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 30 CM - LEFT
267.030R	667.030R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 30 CM - RIGHT
267.032L	667.032L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 32 CM - LEFT
267.032R	667.032R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 32 CM - RIGHT
267.034L	667.034L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 34 CM - LEFT
267.034R	667.034R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 34 CM - RIGHT
267.036L	667.036L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 36 CM - LEFT
267.036R	667.036R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 36 CM - RIGHT
267.038L	667.038L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 38 CM - LEFT
267.038R	667.038R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 38 CM - RIGHT
267.040L	667.040L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 40 CM - LEFT
267.040R	667.040R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 40 CM - RIGHT
267.042L	667.042L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 42 CM - LEFT
267.042R	667.042R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 42 CM - RIGHT
267.044L	667.044L	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 44 CM - LEFT
267.044R	667.044R	LONG PROXIMAL FEMORAL NAILS Ø 09 MM X 44 CM - RIGHT
267.130L	627.130L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 30 CM - LEFT
267.130R	627.130R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 30 CM - RIGHT
267.132L	627.132L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 32 CM - LEFT
267.132R	627.132R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 32 CM - RIGHT
267.134L	627.134L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 34 CM - LEFT
267.134R	627.134R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 34 CM - RIGHT
267.136L	627.136L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 36 CM - LEFT
267.136R	627.136R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 36 CM - RIGHT
267.138L	627.138L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 38 CM - LEFT
267.138R	627.138R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 38 CM - RIGHT
267.140L	627.140L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 40 CM - LEFT
267.140R	627.140R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 40 CM - RIGHT
267.142L	627.142L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 42 CM - LEFT
267.142R	627.142R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 42 CM - RIGHT
267.144L	627.144L	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 44 CM - LEFT
267.144R	627.144R	LONG PROXIMAL FEMORAL NAILS Ø 10 MM X 44 CM - RIGHT
267.230L	667.230L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 30 CM - LEFT
267.230R	667.230R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 30 CM - RIGHT
267.232L	667.232L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 32 CM - LEFT
267.232R	667.232R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 32 CM - RIGHT
267.234L	667.234L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 34 CM - LEFT
267.234R	667.234R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 34 CM - RIGHT
267.236L	667.236L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 36 CM - LEFT
267.236R	667.236R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 36 CM - RIGHT
267.238L	667.238L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 38 CM - LEFT
267.238R	667.238R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 38 CM - RIGHT
267.240L	667.240L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 40 CM - LEFT
267.240R	667.240R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 40 CM - RIGHT
267.242L	667.242L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 42 CM - LEFT
267.242R	667.242R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 42 CM - RIGHT
267.244L	667.244L	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 44 CM - LEFT
267.244R	667.244R	LONG PROXIMAL FEMORAL NAILS Ø 11 MM X 44 CM - RIGHT

Specification

- Available in Stainless steel & Titanium.
- Available in 30cm to 44cm length with 9.0, 10 & 11mm diameter.
- Hip pin 6.2mm, femoral neck 8.0mm, 4.9mm locking bolt is used.
- The nail's distal flexibility and optimal distance between the interlocking bolt and the tip of the nail minimize stress concentration and tension in the femoral shaft.

Indication for use

- The Long femoral nail is used for long stature Indications Low and extended subtrochanteric fractures, Ipsilateral trochanteric fractures, Combination of fractures (trochanteric area / shaft), Pathological fractures Contraindications or combined medial femoral neck fractures.

