



CAT. NO. SS	CAT. NO. TT	LONG PROXIMAL FEMORAL NAILS (135°)
267.667L	-	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 32 CM - LEFT
267.667R	-	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 32 CM - RIGHT
267.668L	667.334L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 34 CM - LEFT
267.668R	667.334R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 34 CM - RIGHT
267.669L	667.336L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 36 CM - LEFT
267.669R	667.336R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 36 CM - RIGHT
267.670L	667.338L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 38 CM - LEFT
267.670R	667.338R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 38 CM - RIGHT
267.671L	667.340L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 40 CM - LEFT
267.671R	667.340R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 40 CM - RIGHT
267.672L	667.342L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 42 CM - LEFT
267.672R	667.342R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 42 CM - RIGHT
267.673L	667.344L	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 44 CM - LEFT
267.673R	667.344R	LONG PROXIMAL FEMORAL NAILS 135° 9 MM X 44 CM - RIGHT
267.682L	667.354L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 34 CM - LEFT
267.682R	667.354R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 34 CM - RIGHT
267.683L	667.356L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 36 CM - LEFT
267.683R	667.356R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 36 CM - RIGHT
267.684L	667.358L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 38 CM - LEFT
267.684R	667.358R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 38 CM - RIGHT
267.685L	667.360L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 40 CM - LEFT
267.685R	667.360R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 40 CM - RIGHT
267.686L	667.362L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 42 CM - LEFT
267.686R	667.362R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 42 CM - RIGHT
267.687L	667.364L	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 44 CM - LEFT
267.687R	667.364R	LONG PROXIMAL FEMORAL NAILS 135° 10 MM X 44 CM - RIGHT
267.698L	667.374L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 34 CM - LEFT
267.698R	667.374R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 34 CM - RIGHT
267.699L	667.376L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 36 CM - LEFT
267.699R	667.376R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 36 CM - RIGHT
267.700L	667.378L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 38 CM - LEFT
267.700R	667.378R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 38 CM - RIGHT
267.701L	667.380L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 40 CM - LEFT
267.701R	667.380R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 40 CM - RIGHT
267.702L	667.382L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 42 CM - LEFT
267.702R	667.382R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 42 CM - RIGHT
267.703L	667.384L	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 44 CM - LEFT
267.703R	667.384R	LONG PROXIMAL FEMORAL NAILS 135° 11 MM X 44 CM - RIGHT
267.706L	667.385L	LONG PROXIMAL FEMORAL NAILS 135° 12 MM X 34 CM - LEFT
267.706R	667.385R	LONG PROXIMAL FEMORAL NAILS 135° 12 MM X 34 CM - RIGHT
267.707L	667.386L	LONG PROXIMAL FEMORAL NAILS 135° 12 MM X 36 CM - LEFT

267.707R	667.386R	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 36 CM - RIGHT
267.708L	667.387L	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 38 CM - LEFT
267.708R	667.387R	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 38 CM - RIGHT
267.709L	667.388L	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 40 CM - LEFT
267.709R	667.388R	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 40 CM - RIGHT
267.710L	667.389L	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 42 CM - LEFT
267.710R	667.389R	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 42 CM - RIGHT
267.711L	667.390L	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 44 CM - LEFT
267.711R	667.390R	LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 44 CM - RIGHT
267.712L		LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 44 CM - LEFT
267.712R		LONG PROXIMAL FEMORAL NAILS 135* 12 MM X 44 CM - RIGHT
267.715L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 32 CM - LEFT
267.715R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 32 CM - RIGHT
267.716L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 34 CM - LEFT
267.716R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 34 CM - RIGHT
267.717L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 36 CM - LEFT
267.717R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 36 CM - RIGHT
267.718L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 38 CM - LEFT
267.718R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 38 CM - RIGHT
267.719L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 40 CM - LEFT
267.719R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 40 CM - RIGHT
267.720L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 42CM - LEFT
267.720R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 42CM - RIGHT
267.721L		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 44CM - LEFT
267.721R		LONG PROXIMAL FEMORAL NAILS 135* 13 MM X 44CM - RIGHT

Specification

- Available in both Stainless steel & Titanium.
- Available in 9, 10, 11, 12 & 13mm diameter from 32 to 44 cm length.
- Compatible with 4.9mm locking bolt, Hip pin 6.2mm, femoral neck screw 8.0mm is used.

Indication for use

The long femoral nail is used for small stature.

- Low and extended subtrochanteric fractures
- Ipsilateral trochanteric fractures
- Combination of fractures (trochanteric area / shaft)
- Pathological fractures