



CAT. NO. SS	CAT. NO. TT	PROXIMAL FEMORAL NAILS 135° (SHORT)
267.658	667.391	PROXIMAL FEMORAL NAILS 135* Ø 09 MM X 24 CM
267.659	667.393	PROXIMAL FEMORAL NAILS 135* Ø 10 MM X 24 CM
267.660	667.395	PROXIMAL FEMORAL NAILS 135* Ø 11 MM X 24 CM
267.661	667.397	PROXIMAL FEMORAL NAILS 135* Ø 12 MM X 24 CM

Specification

- Available in both Stainless steel & Titanium.
- Available in 9, 10, 11 & 12mm diameter with 24 cm length.
- Compatible with 4.9mm locking bolt, Hip pin 6.2mm, femoral neck screw 8.0mm is used.

Indication for use

The short femoral nail is used for small stature.

- Pertrochanteric fractures
- Intertrochanteric fractures
- High subtrochanteric fractures