



CAT. NO. SS	CAT. NO. TT	PROXIMAL FEMORAL NAILS (135°)
267.658	667.391	LONG PROXIMAL FEMORAL NAILS 135* 9 MM X 24 CM
267.659	667.393	LONG PROXIMAL FEMORAL NAILS 135* 10MM X 24 CM
267.660	667.395	LONG PROXIMAL FEMORAL NAILS 135* 11MM X 24 CM
267.661	667.397	LONG PROXIMAL FEMORAL NAILS 135* 12MM X 24 CM

Specification

- Available in both Stainless steel & Titanium.
- Available in 9, 10, 11 & 12mm diameter with 24 cm length.
- Compatible with 4.9mm locking bolt, Hip pin 6.2mm, femoral neck screw 8.0mm is used.

Indication for use

- Proximal Femoral Nails are used to fix pertrochanteric, intertrochanteric, high subtrochanteric and transcervical fractures of the Femur.

