

SHORT PROXIMAL FEMORAL NAILS

130*/135* (SHORT P.F.N.)



CAT. NO. SS	CAT. NO. TT	Short Proximal Femoral Nails 130* /135* (Short P.F.N.)
266.090	666.090	PROXIMAL FEMORAL NAILS Ø 09 MM X 24 CM
266.100	666.100	PROXIMAL FEMORAL NAILS Ø 10 MM X 24 CM
266.110	666.110	PROXIMAL FEMORAL NAILS Ø 11 MM X 24 CM
266.120	666.120	PROXIMAL FEMORAL NAILS Ø 12 MM X 24 CM

Specification

- Available in Stainless steel & Titanium.
- Available in 24cm length with 9.0, 10, 11, 12mm diameter.
- Hip pin 6.2mm, femoral neck screw 8.0mm, 4.9mm locking bolt is used.

Indication for use

- Proximal Femoral Nails are used to fix pertrochanteric, intertrochanteric, high subtrochanteric and transcervical fractures of the Femur.

