

TROCHANTERIC FEMORAL NAILS 130*/135* (TFN)



CAT. NO. SS	CAT. NO. TT	Trochanteric Femoral Nails 130*/135* (TFN)
265.090	669.090	TROCHANTERIC NAILS Ø 09 MM X 18 CM
265.100	669.100	TROCHANTERIC NAILS Ø 10 MM X 18 CM
265.110	669.110	TROCHANTERIC NAILS Ø 11 MM X 18 CM
265.120	669.120	TROCHANTERIC NAILS Ø 12 MM X 18 CM

Specification

- Available in Stainless steel & Titanium.
- Available in 18cm length with 9.0, 10, 11, 12mm diameter.
- Hip pin 6.2mm, femoral neck 8.0mm, 4.9mm locking bolt is used.

Indication for use

- The trochanteric femoral nail (TFN) is a helical blade device which has fins that compact the cancellous bone as it is inserted into bone and may provide improved fracture stabilization characteristics. Implants were placed percutaneously, and compression of the fracture was performed in the majority of fractures.

