



CAT. NO. SS	CAT. NO. TT	TROCHANTERIC NAILS (135°)
267.650	667.401	TROCHANTERIC NAILS 135* Ø 09 MM X 18 CM
267.651	667.403	TROCHANTERIC NAILS 135* Ø 10 MM X 18 CM
267.652	667.405	TROCHANTERIC NAILS 135* Ø 11 MM X 18 CM
267.653	667.407	TROCHANTERIC NAILS 135* Ø 12 MM X 18 CM

Specification

- Available in both Stainless steel & Titanium.
- Available in 9, 10, 11 & 12 mm diameter with 18 cm length.
- Compatible with 4.9mm locking bolt, Hip pin 6.2mm, femoral neck screw 8.0mm.

Indication for use

- Treating Fractures of the Proximal Femur including Pertrochanteric Fractures, Intertrochanteric Fractures, Basal Neck Fractures.