



CAT. NO.	METAL	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW
692.742	TT	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW 12.0MM
692.744	TT	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW 14.0MM
692.746	TT	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW 16.0MM
692.748	TT	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW 18.0MM
692.750	TT	3.0MM TITANIUM CERVICAL SPINE LOCKING SCREW 20.0MM

Specification

- Available in Titanium.
- Available in 12,14,16,18 & 20mm size with 3.0mm Diameter.

Indication for use

- It is used for the internal fixation in the treatment of instabilities associated with fractures/dislocations, degenerative diseases, tumours and partial or total spondylectomy.