



CAT. NO.	METAL	LCP Variable Angle Cervical Screw 3.5mm
690.312	TT	LCP Variable Angle Cervical Screw 3.5mm x 12.0mm
690.314	TT	LCP Variable Angle Cervical Screw 3.5mm x 14.0mm
690.316	TT	LCP Variable Angle Cervical Screw 3.5mm x 16.0mm
690.318	TT	LCP Variable Angle Cervical Screw 3.5mm x 18.0mm
690.320	TT	LCP Variable Angle Cervical Screw 3.5mm x 20.0mm

## Specification

- LCP Variable Angle Cervical Screw is available in Titanium.
- Available in 12,14,16,18 & 20mm size with 3.5mm Dia.

## Indication for use

- It provides the surgeon some flexibility on the trajectory of the locking screws in order to facilitate fixation of variable fracture patterns and anatomy.