



CAT. NO.	METAL	CERVICAL PLATE WITH INBUILT CAGE
692.605	TT	Cervical Plate with Inbuilt Cage Ø5.0mm
692.606	TT	Cervical Plate with Inbuilt Cage Ø6.0mm
692.607	TT	Cervical Plate with Inbuilt Cage Ø7.0mm
692.608	TT	Cervical Plate with Inbuilt Cage Ø8.0mm
692.609	TT	Cervical Plate with Inbuilt Cage Ø9.0mm
692.610	TT	Cervical Plate with Inbuilt Cage Ø10.0mm

Specification

- Cervical Plate with Inbuilt Cage available in Titanium.
- Available in 5,6,7,8,9 & 10mm Diameter.

Indication for use

- It can increase fusion rates, maintain or improve cervical sagittal alignment and primary and secondary stability.