

INNILOCK MONOAXIAL SCREW 5.5MM (FOR 6MM ROD)



CAT. NO.	METAL	INNILOCK Monoaxial Screw 5.5mm (For 6mm rod)
643.130	TT	INNILOCK Monoaxial Screw- Ø 5.5 mm x 30 mm
643.135	TT	INNILOCK Monoaxial Screw- Ø 5.5 mm x 35 mm
643.140	TT	INNILOCK Monoaxial Screw- Ø 5.5 mm x 40 mm
643.145	TT	INNILOCK Monoaxial Screw- Ø 5.5 mm x 45 mm
643.150	TT	INNILOCK Monoaxial Screw- Ø 5.5 mm x 50 mm

Specification

- Available length 30,35,40,45 & 50mm with 5.5mm dia.

Indication for use

- The 5.5 MM Monoaxial screw is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies).
- Pedicle screw fixation is limited to skeletally mature patients.
- The 5.5 MM Monoaxial Screw is a non-cervical spinal fixation device intended for use as a posterior pedicle screw.