

INNILOCK MONOAXIAL SCREW 6.5MM (FOR 6MM ROD)



| CAT. NO. | METAL | INNILOCK Monoaxial Screw 6.5mm (For 6mm rod) |
|----------|-------|---|
| 643.230 | TT | INNILOCK Monoaxial Screw - Ø 6.5 mm x 30 mm |
| 643.235 | TT | INNILOCK Monoaxial Screw - Ø 6.5 mm x 35 mm |
| 643.240 | TT | INNILOCK Monoaxial Screw - Ø 6.5 mm x 40 mm |
| 643.245 | TT | INNILOCK Monoaxial Screw - Ø 6.5 mm x 45 mm |
| 643.250 | TT | INNILOCK Monoaxial Screw - Ø 6.5 mm x 50 mm |

Specification

- Available length 30,35,40,45 & 50mm with 6.5mm dia.

Indication for use

- The 6.5 MM Monoaxial screw is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies).
- Pedicle screw fixation is limited to skeletally mature patients.
- The 6.5 MM Monoaxial Screw is a non-cervical spinal fixation device intended for use as a posterior pedicle screw.

